

## 40-700 ESCR-Tester

#### **Standards**

**ASTM D 1693** 







10-Station-Device

50-Station-Device

Tools

### **Application**

The ESCR test measures the susceptibility of ethylene plastics to environmental stress-cracking. Under certain conditions of stress and in the presence of environments such as soaps, wetting agents, oils or detergents, ethylene plastics may exhibit mechanical failure by cracking. Bent specimens of the plastic, each having a controlled imperfection on one surface, are exposed to the action of a surface-active agent. The proportion of the total number of specimens that crack in a given time is observed.

#### **Features**

The ESCR tester comprises a temperature bath with a test tube racks holding up to 10/50 test tubes. Each test tube has space for 10 specimens. For sample preparation, a notching unit with natching blade is supplied. The notched specimens are pre-stressed using the bending clamp and inserted into the specimen holder using the transfer tool.

Up to 10 test tubes are placed in a rack, which can be removed from the ESCR test bath for sample preparation and inserted into the bath for testing. The 50 Station unit can hold up to 5 sample racks and can thus be configured as a 10/20/30/40/50 station unit. 30 samples can be viewed through the front window and up to 50 can be tested simultaneously.

A unique feature of the Coesfeld ESCR tester is the large glass front with integrated LED light to inspect samples while sample remain in the bath. This makes measurement time efficient and clean.

#### **Technical Data**

	10-Station-Device	50-Station-Device
Dimensions (WxDxH)	526 x 440 x 476 mm	1331 x 440 x 476 mm
Bath Volume	~24	~72
Weight	~35 kg (empty)	~75 kg (empty)
Temperature range	RT + 10 °C to 150 °C	
Temperature stability	± 0.1 °C	
Mains / Power	230 V, 50 Hz (alternative on request) / 2000 W	



# List of items

Item	Description		
40-700-011		10 Station ESCR Tester (up to 10 specimen)	
40-700-010		50 Station ESCR Tester (up to 50 specimen)	
9-153-123		Lid, covering bath when sample rack is not inserted/ 1 lid per test rack	
9-953-010		Test Rack for holding up to 10 test glasses/ min. 1 rack required per unit/ max.5 racks possible for 50 Station tester	
40-700-013		Test glass, outer $\emptyset$ 36 mm, inner $\emptyset$ 32 mm, length = 225 mm, 1 per sample holder	
40-700-002	35	Test tube tong to hold test glass / 1 per unit	
40-700-003		Cork stopper, 1 per test glass	
40-700-004	See also de la companya de la compan	Brass specimen holder for 10 samples / 1 per test glass	
9-953-004		Nicking jig to cut sample / 1 per unit	
40-700-006		Gauge to measure notch depth / 1 per nicking jig	
40-700-005		Notching Blades for nicking jig (5/pack) (1 Blade can be used for notching approx. 100 samples)	
40-700-012		Notching Blades for nicking jig (200/pack) (1 Blade can be used for notching approx. 100 samples)	
9-953-005		Bending clamp to simultaneously bend 10 sample / 1 per unit	
9-953-006		Transfer tool to transfer sample from bending clamp to brass holder / 1 per unit	